

2001 Scrubber Unit

Parts List


1000 8925 / Rev. 1.3



Note: Screws, nuts and washers are shown but not designated in this part list!

Designation: 2001 Scrubber Unit
Drawing: BLOC068e
Effective from serial no: _____ to serial no: _____

Pos	Part number	Designation	Remarks
1	1001 2448	Hose connection	
2	1000 6757	O-rings set	
3	1000 6707	Switch	
4	1565 0233	Power inlet	
5	1564 0080	Fuse holder	
6	1564 0123	Fuse holder	
7	1001 2395	Scrubber pump cpl.	
8	1001 2296	Magnetically actuated valve	
9	1584 0002	Ball valve	
10	1565 0116	Contact housing	
11	1586 0009	Hose adaptor	
12	1562 0022	Solid state relay	
13	1568 0029	Transformer	
14	1544 0017	Vibration absorber	
15	1582 0162	Tubing (viton)	
16	1564 0097	Fuse	
17	1000 6390	Lid	
18	1000 6401	Flask	



BUCCO®
1000 6015

SCRUBBER

This page is intentionally not used.